



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare Gynaecology Couches (SM8563) (New)

- New,
- Made in UK,
- Electric control of height and tilt angle,
- Adjustable angle of backrest and leg support controlled manually with gas spring,
- Compact frame design with retractable wheel system (RWD),
- Crescent shaped seat design,
- Basics:
 - Paper Roll holder (6005),
 - Regulated patient leg support,
 - Debris tray,
- Base frame powder coated with antibacterial layer,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 147 cm (without leg support),
 - Width: 65 cm,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from horizontal to 80 °,
 - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
 - Removable headrest - pillow (6007),
 - Fully adjustable patient leg supports (6043),
 - Extra folding patient leg support (6051),
 - Short base cover (6052),
 - Battery back-up (6057),
 - Standard operators chair (6100),
 - Foot switch (6200),
- Technical Passport,
- 24-month warranty for electric



Model 6100



Model 6051



Model 6043

