



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medicare Gynaecology Couches (SM8573) (New)**

- New,
- Made in UK,
- Electric control of height and tilt angle,
- Supplied with matching operators chair and head cushion,
- Compact frame design with retractable wheel system (RWD),
- Crescent shaped seat design,
- Basics:
  - Paper Roll holder (6005),
  - Regulated patient leg support,
  - Debris tray,
- Base frame powder coated with antibacterial layer,
- TiMotion electric linear actuators,
- Technical Specifications:
  - Length: 147 cm (without leg support),
  - Width: 65 cm,
  - Safe working load and lifting capacity of 240 kg (528lbs),
  - Adjustable angle electric backrest from horizontal to 80 °,
  - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
  - Removable headrest - pillow (6007),
  - Fully adjustable patient leg supports (6043),
  - Extra folding patient leg support (6051),
  - Short base cover (6052),
  - Battery back-up (6057),
  - Standard operators chair (6100),
  - Triple foot switch (6209),
- Technical Passport,
- 24-month warranty for electric



Model 6100



Model 6051



Model 6043



Height range 52cm to 106cm

