



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medicare GP Gynaecology Couches (SM8543) (New)

- New,
- Made in UK,
- Electric height, tilt and backrest with manual head section,
- Retractable wheel system,
- Vertical lifting design with stability throughout,
- Crescent shaped seat design,
- Basic equipment:
  - Breathing hole and plug for prone treatments (6004),
  - Paper Roll handle (6005),
  - Lithotomy stirrups (6045),
  - Debris tray,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
  - Length: 196 cm,
  - Width: 70 cm,
  - Safe working load and lifting capacity of 240 kg (528 lbs),
  - Adjustable angle electric backrest from 0 ° to 80 °,
  - Power assisted adjustable head section from -25 ° to 50 °,
  - Excellent height range from 50 cm to 103 cm,
- **Additional Equipment (extra paid):**
  - Removable headrest (6007),
  - Lithotomy stirrups (6045),
  - Protective Base Shrouding (6053),
  - Battery Back-up (6057),
  - Standard operators chair (6100),
  - Triple foot switch (6209),
- Technical Passport,
- 24-month warranty for electric



Model 6045



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATEK)

