



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare GP Gynaecology Couches (SM8543D) (New)

- New,
- Made in UK,
- Deluxe Model complete with base cover, matching operators chair and head support cushion,
- Retractable wheel system
- Vertical lifting design with stability throughout,
- Crescent shaped seat design,
- Basic equipment:
 - Breathing hole and plug for prone treatments (6004),
 - Paper Roll handle (6005),
 - Lithotomy stirrups (6045),
 - Debris tray,
 - Protective Base Shrouding (6053),
 - Removable headrest (6007),
 - Standard operators chair (6100),
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 196 cm,
 - Width: 70 cm,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0 ° to 80 °,
 - Power assisted adjustable head section from -25 ° to 50 °,
 - Excellent height range from 50 cm to 103 cm,
- **Additional Equipment (extra paid):**
 - Lithotomy stirrups (6045),
 - Protective Base Shrouding (6053),
 - Battery Back-up (6057),
 - Triple foot switch (6209),
- Technical Passport,
- 24-month warranty for electric




Model 6045



TAPICERKI STANDARD +[ZA DOPŁATK.]



Height range 50cm to 103cm



TAPICERKI STANDARD [W CENIE]



