



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl




Examination Table Therapy 3 Section (ST3568) (New)

- New,
- Made in UK,
- Electric height elevation,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Hand switch fitted as standard on electric models,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 188 cm,
 - Width: 65 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
 - Adjustable headrest: -25° to 75°,
 - Adjustable Backrest: 0° to 80°,
 - Excellent height range from 45 cm (18") to 98 cm (38"),
- **Additional Equipment (extra paid):**
 - Perimeter foot switch (EPFS) (ST6001),
 - Side Support Rails (ST6002),
 - Side Support Cushions (ST6002P),
 - Paper towel holder (ST6005),
 - Upper limb support (ST6006),
 - Removable head/knee support (ST6007),
 - Matching therapy chair or stool,
 - Single foot switch (ST6200),
 - Double foot switch (ST6239),
- Technical Passport,
- 24-month warranty for electrics and gas springs



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATEK]

