



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Table Therapy 3 Section Drop (ST3566) (New)

- New,
- Made in UK,
- Electric Height Elevation,
- The power assisted backrest and foot sections incorporate anti-trap design gas struts,
- Head section standard with breathing hole and plug,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Hand switch fitted as standard on electric models,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 191 cm,
 - Width: 65 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37.7 stone),
 - Adjustable angle backrest from horizontal to +80 degrees positive,
 - Adjustable angle head section from +10 to -90 degrees negative,
 - Excellent height range from 48 - 101 cm (19" - 39"),
- **Additional Equipment (extra paid):**
 - Perimeter foot switch (EPFS) (ST6001),
 - Side Support Rails (ST6002),
 - Side Support Cushions (ST6002P),
 - Paper towel holder (ST6005),
 - Upper limb support (ST6006),
 - Removable head/knee support (ST6007),
 - Matching therapy chair or stool,
 - Single foot switch (ST6200),
 - Double foot switch (ST6239),
- Technical Passport,
- 24-month warranty for electrics and gas springs



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATKĄ]