



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)


Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Examination Table Deluxe Therapy Non - Drainage (ST3349) (New)

- New,
- Made in UK,
- Electric height elevation via integral perimeter foot switch,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- Deluxe Therapy Couch with head design , Perimeter foot switch operation for electric models (ST6001) - armrests,
- TiMotion electric linear actuators,
- Technical Data:
  - Length: 195 cm,
  - Width: 65 cm,
  - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
  - Adjustable headrest: -25° to 40°,
  - Adjustable leg plates: 0° to 90°,
  - Excellent height range from 46 - 99 cm (19"-39"),
- **Additional Equipment (extra paid):**
  - Paper towel holder (ST6005),
  - Removable head/knee support (ST6007),
  - Protective Base Shroud (ST6053),
  - Matching therapy chair or stool,
- Technical Passport,
- 24-month warranty for electrics and gas springs



TAPICERKI STANDARD [W CENIE ]



TAPICERKI STANDARD +[ZA DOPŁAT.]

