



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Table Therapy Hygiene Table (ST1661) (New)

- New,
- Made in UK,
- Hygiene Therapy Table - electric, large mobile wheel design,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 184 cm,
 - Width: 65 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
 - Excellent height range from 43 cm (17") to 97 cm (37"),
- **Additional Equipment (extra paid):**
 - Side support rails (ST6002),
 - Side support cushions (ST6002P),
 - Paper towel holder (ST6005),
 - Stretch towelling covers (ST6010), (ST6011),
 - Removable head/knee support (ST6007),
 - Protective base shrouding (ST6053),
 - Single foot switch (ST6200),
 - Double foot switch (ST6239),
 - Matching therapy chair or stool,
- Technical Passport,
- 24-month warranty for electrics and gas springs

