



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)


Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Examination Table Sterling Changing Short Table (SX1057) (New)

- New,
- Made in UK,
- Electric Elevation,
- 100 mm diameter individual locking castors,
- Basic equipment:
  - Side support rails,
  - Remote control,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- Technical Specifications:
  - Length: 184 cm,
  - Width: 63 cm,
  - Safe working load and lifting capacity (from minimum height) of 200 kg (440 lbs = 31.5 stone),
  - Excellent height range from 46 cm to 99 cm,
- **Additional Equipment (extra paid):**
  - Side support rails cushions (Pair) (6003P),
  - Paper Roll holder (6005),
  - Foot switch (6200),
- Technical Passport,
- 24-month warranty for electrics



Width  
63cm



184cm  
Model SX1057



TAPICERKI STANDARD [W CENIE ]



TAPICERKI STANDARD +[ZA DOPŁATEK]

