



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Table Sterling Changing Short Table (SX1057) (New)

- New,
- Made in UK,
- Electric Elevation,
- 100 mm diameter individual locking castors,
- Basic equipment:
 - Side support rails,
 - Remote control,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 184 cm,
 - Width: 63 cm,
 - Safe working load and lifting capacity (from minimum height) of 200 kg (440 lbs = 31.5 stone),
 - Excellent height range from 46 cm to 99 cm,
- **Additional Equipment (extra paid):**
 - Side support rails cushions (Pair) (6003P),
 - Paper Roll holder (6005),
 - Foot switch (6200),
- Technical Passport,
- 24-month warranty for electrics



Width
63cm



184cm
Model SX1057



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATEK]

