



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Chair Surgeons & Sonographers (MC6158) (New)

- New
- Made in UK
- Surgeons chair - standard model
- Height range 45cm to 59cm.
- Wide choice of models available
- Ergonomically designed seat
- Compact backrest for improved scanning procedures
- Seat and backrest upholstered on both sides
- Weight capacity and safe working load of 150Kg
- Anti-bacterial and waterproof vinyl upholstery
- Restraint castors fitted as standard
- **Dostępne akcesoria (za dopłatą):**
  - A set of anti-static castors (6178),
  - Load loading castors automatically locked (6455),
  - Manual Locking castors (6172),
  - Glides (static feet) (6173),
  - Drop down arm supports (6175),
  - Foot Height Control (FHC) (6174),
  - Foot support ring (6176),
- 24 months warranty



TAPICERKI STANDARD + (ZA DOPŁATĄ)



TAPICERKI STANDARD (W CENIE)

