



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Phlebotomy Chair (MC6154) (New)

- New
- Made in UK
- The chair is a compact and space saving alternative to electrically operated chairs and couches.
- Supplied with a sturdy base on glides (static feet) or with several castor options.
- Integral arm supports are fitted as standard to assist the patient in mounting and dismounting the chair.
- The chair has a height range of 49cm to 63cm
- Max load: 200kg,
- 24 months warranty

