



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Medical Stool Coburg Dentalift 22015 (New)

- New
- Made in Germany
- dental operator's stool (according to DIN EN ISO 7493)
- base made of stainless steel with 5 twin castors for hard floors
- chrome-plated foot ring
- height adjustment using a safety gas spring by means of an overlying foot-controlled disc
- seat and back rest anatomically-shaped and easily adjustable by means of a lever
- comfortable upholstery as per our sample cards (special colours subject to a surcharge)
- Seat Height:
 - L: 42 - 56 cm,
 - T: 48 - 68 cm,
- **Additional Equipment (paid):**
 - metal parts in WHITE LINE - RAL 9002, powder-coated
 - operator's armrest 2014, swivel and height adjustable, mountable at right or/and left side
 - special heights on request
- Technical Passport
- 24 months warranty

