



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medical Stool Coburg Medicalift 3014 (New)**

- New
- Made in Germany,
- Special surgical stool suitable for all areas of surgery, particularly microsurgery (in accordance with DIN EN ISO 7493)
- Stainless steel safety base with 5 anti-static twin rollers
- Chrome foot ring
- Upper part can be rotated and locked into position
- Seat and back rest anatomically shaped and easily adjustable by means of a lever
- Comfortable upholstery, anti-static, black
- Height adjustment via hydraulic pump lever; lift height: 19 cm
- Seat Height: 53 - 72 cm
- **Additional Equipment (paid):**
  - central stop bar
  - operator armrests 2016
    - 1 pair of flat stainless steel armrests, can be sterilized
    - adjustable height and can be rotated on all sides by means of a ball joint
    - available with washable upholstery if required
- Technical Passport,
- 24 months warranty

