



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Medical Stool Coburg Medicalift 3030 (New)

- New
- Made in Germany,
- Special surgical chair suitable for all areas of surgery, particularly microsurgery
- Stainless steel safety base with 5 anti-static twin rollers
- Chrome foot ring
- Upper part can be rotated and locked into position
- Back rest can be adjusted by 7cm
- Comfortable upholstery, anti-static, black
- Height adjustment with maintenance-free hydraulic pump lever, lifting height: 19 cm
- Seat height: 53 - 72 cm
- **Additional Equipment (paid):**
 - pair of upholstered operator armrests 2030
 - central stop bar
- Technical Passport,
- 24 months warranty



	<p>KOŁOWY PRĄT Wzrost regulowany i bezpieczny wspinaczki. Idealny dla mikrochirurgii ze sterowaniem przez pokrętkę hydrauliczną. Nie posiada podłokietników.</p> <p>Modely ze sterowaniem ręcznym: 3030, 3030A, 3030A, 3030A, 3030A, 3030A, 3030A</p>
	<p>KOŁOWY PRĄT Pudełkowy ze sterowaniem i regulacją (opcjonalnie). Wyposażony w amortyzatory regulowane. Dostępny w wersji bez podłokietników.</p> <p>Modely ze sterowaniem ręcznym: 3030, 3030A, 3030A, 3030A, 3030A, 3030A</p>
	<p>KOŁOWY PRĄT Wzrost regulowany i bezpieczny wspinaczki. Wyposażony w amortyzatory regulowane. Dostępny w wersji bez podłokietników.</p> <p>Modely ze sterowaniem ręcznym: 3030, 3030A, 3030A, 3030A, 3030A, 3030A</p>