



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medical Stool Coburg Medicalift 3030 (New)

- New
- Made in Germany,
- Special surgical chair suitable for all areas of surgery, particularly microsurgery
- Stainless steel safety base with 5 anti-static twin rollers
- Chrome foot ring
- Upper part can be rotated and locked into position
- Back rest can be adjusted by 7cm
- Comfortable upholstery, anti-static, black
- Height adjustment with maintenance-free hydraulic pump lever, lifting height: 19 cm
- Seat height: 53 - 72 cm
- **Additional Equipment (paid):**
  - pair of upholstered operator armrests 2030
  - central stop bar
- Technical Passport,
- 24 months warranty

