



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**Medical Stool Coburg Medicalift 22011 (New)**

- New
- Made in Germany,
- Multi-purpose stool suitable for surgery and outpatient departments
- Safety base made of stainless steel with 5 anti-static twin castors for hard floors
- Chrome-plated foot ring
- Height adjustment using a safety gas spring controlled by means of stainless steel foot-controlled disc
- Seat height: 60 - 83 cm
- Standard version: black saddle seat
- **Additional Equipment (paid):**
  - backrest mod. 2\*
  - rounded seat Ø35cm\*
  - rounded seat Ø38cm\*
  - moulded seat 1103\*
  - comfortable upholstery as per our sample card
  - special colours subject to a surcharge
- Technical Passport,
- 24 months warranty

