



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Medical Stool Coburg Medicalift 22011 (New)

- New
- Made in Germany,
- Multi-purpose stool suitable for surgery and outpatient departments
- Safety base made of stainless steel with 5 anti-static twin castors for hard floors
- Chrome-plated foot ring
- Height adjustment using a safety gas spring controlled by means of stainless steel foot-controlled disc
- Seat height: 60 - 83 cm
- Standard version: black saddle seat
- **Additional Equipment (paid):**
 - backrest mod. 2*
 - rounded seat Ø35cm*
 - rounded seat Ø38cm*
 - moulded seat 1103*
 - comfortable upholstery as per our sample card
 - special colours subject to a surcharge
- Technical Passport,
- 24 months warranty

