



## Medical Stool Coburg Medicalift 22012 (New)

- New
- Made in Germany,
- Multi-purpose stool suitable for surgery, anaesthesia and outpatient departments
- Stainless steel security base with 5 twin castors for hard floors
- Height adjustment using a safety gas spring controlled by means of a stainless steel foot-controlled disc
- Rounded seat with anti-static black upholstery
- Seat Height:
  - H200: 47 - 66cm
  - H250: 52 - 76cm
- **Additional Equipment (paid):**
  - back rest with anti-static black upholstery
  - special colours subject to a surcharge
  - central circular stop bar operated by placing your feet on the footplate
- Technical Passport,
- 24 months warranty

