



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Coburg Exaflex 6126 operating table (New)

- A brand new product,
- Made in Germany,
- **The basic version of the table includes:**
 - A headrest (standard),
 - A footrest,
 - Driving wheels system,
- Table used for examination in maxillofacial surgery, ENT, ophthalmology, gynecology, dermatology, urology, microsurgery and other medical departments,
- Stability of the table is provided by a spacious base plate (700 x 500 mm) with a stainless steel cover,
- Table top (2-part), backrest and footrest made of cast aluminum,
- Handrails on both sides: seat, backrest and footrest made of stainless steel,
- Easily pull-out headrest,
- Metal parts of the table powder-coated RAL 9007,
- Regulation options:
 - Height regulation 60 - 90 cm by electro-mechanical lever, controlled by foot pedal,
 - Vertical tilt Trendelenburg position (20°) by the use of gas spring,
 - Backrest regulation + 90 ° to -10 ° by the use of gas spring,
 - Footrest regulation + 4 ° C to -90 ° by the use of gas spring,
 - Complete footrest with removable seat extension by loosening 2 knobs (gynecological model),
 - Headrest extendable by 20 cm (clamping mechanism),
- Gynecological version:
 - Without footrest and seat extension (discounted from table price),
 - Drip tray,
 - Leg holders with clamps,
- **Additional options (for extra fee):**
 - Antistatic upholstery,
 - Double thickness of upholstery,
 - Expansion kit up to 55 cm, including upholstery,
- **Additional accessories to the table, visible in the photo (for extra fee),**
 - **24 months for the international market,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **24 months for the domestic market (Poland),**
- Financing possibilities (Poland only): Installments, leasing, loan,

