



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Coburg Exaflex 3000 operating table (New)**

- A brand new product,
- Made in Germany,
- **The basic version of the table has:**
  - A headrest (standard),
  - A footrest (pojedynczy),
- Portable multifunctional surgical table for general use in hospitals and ambulatory care units,
- Multi-purpose versatility thanks to optimal patient positioning and 5,6, or 8-fold tabletop segments.
- The table comes standard with chrome-nickel steel rails located on both sides of the seat, backrest and footrest,
- Upholstery lining with a thickness of 40 mm made of electrically conductive material in black (Other colors and upholstery made of non-conductive material also available),
- Height regulation by a hydraulic pump lever,
- Regulation of backrest, Trendelenburg and reverse-Trendelenburg positions by the use of gas spring,
- Lateral tilting of the table by a crank mechanism,
- Movable base made of stainless steel with large dual rollers with central locking mechanism,
- Column cover made of chrome-nickel steel,
- The table has the possibility of mounting an X-ray cassette,
- Maximum load capacity of 160kg,
- Metal parts of the table powder coated RAL 9007,
- **Additional options (for extra fee):**
  - Sliding table top,
  - Footrest (splittable)
- **Additional accessories for the table, visible in the photo (for extra fee),**
  - **24 months for the international market,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
  - **24 months for the domestic market (Poland),**
- Financing possibilities (Poland only): Installments, leasing, loan,

