



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LED NGP 11 SLIM negatoscope (New)

- Brand new product,
- Made in Poland,
- The LED NGP SLIM series negatoscopes are perfect for the use in operating rooms as well as in treatment rooms, offices specializing in mammography diagnostics, and in many other healthcare facilities,
- The NGP SLIM series is an offer of negatoscopes which use LED technology,
- An attractive design of these devices causes significant increase in the aesthetic value of the arrangement of each surgery,
- Among a wide range of products, models designed to be placed on a desk or for wall mounting can be chosen,
- Advantages of SLIM series negatoscopes:
 - The undoubted advantage of this negatoscope is the original case (modern style),
 - Ultra thin - only 35 mm,
 - The function of automatic switching on and off during photographic film change ensures that the radiologist's eyes do not come into contact with the dazzling light,
 - Light source (LED technology) - works more than 50,000 hours,
 - The negatoscope emits light with an intensity of 6000 cd/m² - 19,000.00 lux regulated in the range between 10-100% of the maximum value,
 - The operation of the device does not cause light flickering effect,
 - The front part of the case is made of acrylic with a thickness of 3 mm, which ensures consistent light diffusion (at > 95%),
 - The highest quality of the device is confirmed by certificates,
 - Power consumption is at the level of 45% of a standard negatoscope (energy-saving model),
- Mounting types of the negatoscope:
 - Wall-mounted,
 - Desk-based,
- Technical data:
 - Power supply: 230 V, 50; 60 Hz,
 - Power consumption: 65 W,
 - Light intensity: 6000 cd/m² - 19,000.00 lux,
 - Light intensity regulation: Yes,
 - Screen dimensions: 36 x 43 cm,
 - Automatic light switch off: in the handle,
 - The number of frames: 1,
 - Illumination uniformity: > 95%, > 95%,
 - Mounting: wall, desk base,
 - Dimensions: 430 x 520 x 35 mm (ultra thin),
 - Weight: 4,5 kg,
 - LED technology,
- It has a valid Technical Passport,
- Guarantee: 12 months,

