



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **NGP 800 P vertical negatoscope (New)**

- A brand new product,
- Made in Poland,
- The negatoscope provides an excellent X-ray image that, thanks to its precision, shows even the smallest lesions in soft tissues and any micro degeneration of other kinds,
- **Advantages of using the vertical NGP 800 P negatoscope:**
  - Uniform distribution of light across the screen,
  - Even and uniform distribution of daylight (Day Light) emission,
  - Easy and fast service, including replacement of light sources,
  - Highly regulated light intensity,
  - Durable case,
- **Technical data:**
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 550 VA,
  - Screen luminance: 4100 cd/m<sup>2</sup>±25%,
  - Luminance regulation: optional,
  - Screen dimensions: 72 x 117 cm,
  - Number of frames: 1,
  - Protection class: I,
  - Dimensions: 1270 x 760 x 112 mm,
  - Weight: 35 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

