



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 10 hf negatoscope (New)

- A brand new product,
- Made in Poland,
- The NGP 10 Negatoscope is a product that guarantees high uniformity of light emission over the entire screen surface,
- It has a uniform screen and easy serviceability,
- What's more, it has an ever-reliable roller mount and a durable attractive powder-coated case,
- It has successfully obtained CE declaration along with registration in the Register of Medical Devices and Biocidal Products,
- **Advantages of NGP 10 hf negatoscope:**
 - The negatoscope emits day light (Day Light) with a color temperature of 6500K,
 - High uniformity of light emission across the screen (4xTLD15W),
 - Uniform screen,
 - Easy maintenance (replacement of light sources),
 - Durable attractive powder coated case,
 - All required certificates including: entry in the Register of Medicinal Products, Medical Devices and Biocidal Products CE Declaration of Conformity ISO 13485
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 65 W,
 - Screen luminance: 6000 cd/m²±25%,
 - Luminance regulation: No,
 - Screen dimensions: 36 x 43 cm,
 - Number of frames: 1,
 - Protection class: I,
 - Dimensions: 380 x 510 x 92 mm,
 - Weight: 4,9 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

