



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 20 hf negatoscope (New)

- A brand new product,
- Made in Poland,
- The NGP 20 HF has multiple certifications that prove its superior quality,
- It has a smooth luminance with wide range from 100 up to 50%,
- **Advantages of NGP 20 hf negatoscope:**
 - This model emits daylight with a color temperature of 6500K,
 - Perfect light distribution over the entire screen surface,
 - Uniform screen,
 - Easy installation and maintenance,
 - Durable and modern powder coated case,
 - This device has many certificates that prove its top quality,
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 130 W,
 - Screen luminance: 6000 cd/m²±25%,
 - Luminance regulation: No,
 - Screen dimensions: 72 x 43 cm,
 - Number of frames: 2,
 - Protection class: I,
 - Dimensions: 740 x 510 x 92 mm,
 - Weight: 9,5 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

