



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **NGP 40 hf negatoscope (New)**

- A brand new product,
- Made in Poland,
- The NGP 40 HF model is an ideal product for use in hospitals, clinics, outpatient clinics, treatment rooms, emergency rooms, and many other healthcare facilities,
- **Advantages of NGP 40 hf negatoscope:**
  - Four frames for viewing X-ray negatives,
  - Quick and easy installation,
  - Modern case design,
  - Color temperature of emitted daylight 6500K,
  - Case type: powder coated.
  - Numerous certificates,
  - Excellent screen illumination uniformity parameters
- **Technical data:**
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 260 W,
  - Screen luminance: 6000 cd/m<sup>2</sup>±25%,
  - Luminance regulation: No,
  - Screen dimensions: 144 x 43 cm,
  - Number of frames: 4,
  - Protection class: I,
  - Dimensions: 1460 x 510 x 92 mm,
  - Weight: 18,7 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

