



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 40 hf negatoscope (New)

- A brand new product,
- Made in Poland,
- The NGP 40 HF model is an ideal product for use in hospitals, clinics, outpatient clinics, treatment rooms, emergency rooms, and many other healthcare facilities,
- **Advantages of NGP 40 hf negatoscope:**
 - Four frames for viewing X-ray negatives,
 - Quick and easy installation,
 - Modern case design,
 - Color temperature of emitted daylight 6500K,
 - Case type: powder coated.
 - Numerous certificates,
 - Excellent screen illumination uniformity parameters
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 260 W,
 - Screen luminance: 6000 cd/m²±25%,
 - Luminance regulation: No,
 - Screen dimensions: 144 x 43 cm,
 - Number of frames: 4,
 - Protection class: I,
 - Dimensions: 1460 x 510 x 92 mm,
 - Weight: 18,7 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

