



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NGP 80 hf negatoscope (New)**

- A brand new product,
- Made in Poland,
- The NGP 80 HF eight-frame model is a device that allows for accurate imaging and the ability to see even the smallest changes in soft tissues and any micro degeneration of other types, which can be crucial in the later diagnosis of cancerous conditions,
- **Advantages of NGP 80 hf negatoscope:**
  - All required certifications,
  - Modern case design,
  - Eight-frame device,
  - Emit daylight "Day Light",
  - Have easy installation and service,
  - Have a durable and modern powder-coated case,
  - Have evenly illuminated screen,
- **Technical data:**
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 520 W,
  - Screen luminance: 6000 cd/m<sup>2</sup>±25%,
  - Luminance regulation: No,
  - Screen dimensions: 144 x 96 cm,
  - Number of frames: 8,
  - Protection class: I,
  - Dimensions: 1460 x 960 x 92 mm,
  - Weight: 36 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

