



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



X-Ray View Box NGP 21 m Z (New)

- New,
- Made in Poland,
- edge sharp collimation adjustment even of small image areas,
- brilliant daylight, approx. 6500 K,
- X-ray film viewer with high frequency converters high frequency,
- very quick ignition of fluorescent lamps,
- no flickering effect,
- step-less luminance adjustment within the range from 10% to 100%,
- energy saving up to 25%,
- collimation adjustment combined with luminous density control significantly improves identification of even the tiniest details (microcalcifications), a feature of decisive diagnostic relevance in mammography,
- excellent uniformity of the luminous density,
- brilliant daylight, approx. 6500 K,
- safe X-ray film holder on the whole round of the screen,
- neatly styled sheet steel housing with screwless panel, RAL 9002, structured powder-coated finish, resistant to commercially available disinfectants,
- mobile film holder enabling viewing of standard radiograms,
- X-ray film viewers have all necessary certificates,
- Technical Data:
 - power supply: 230 V 50 Hz,
 - power requirement: 260 VA,
 - luminance: >7000 cd / m²,
 - luminance adjustment: YES,
 - range of adjustment: 10 - 100%,
 - screen size: 48 x 36 cm,
 - frame quantity: 1,
 - class of protection: I,
 - dimensions: 680 x 720 x 120 mm,
 - mass: 55 kg,
- Technical Passport,
- 12-month warranty

