



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Sliding magnifier on the arm (New)

- A brand new product,
- Made in Poland,
- The movable magnifier system can be used in any loupe negatoscope for mammography. It features a precise lamp movement mechanism,
- Technical data;
 - Diameter: 11 cm,
 - Magnification: 2x, 4x,
- Guarantee 12 months,

