



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NBVE 110/55 SL germicidal lamp (New)**

- A brand new product,
- Made in Poland,
- It has a PZH certificate (National Institute of Hygiene),
- NBVE 110/55 SL dual-purpose flow germicidal lamp with two operating time meters,
- The use of energy-saving flow bactericidal and virucidal lamps of the NBVE series is one of the most effective methods of disinfection,
- This device can be used in two ways: when patients and staff are inside (then the airflow function of the lamp is used, and thus the biocidal UV-C radiation does not escape) or when no one is in the room,
- In this case, external irradiators are activated, which emit radiation with a wavelength of 253 nm into the environment, and thus effectively inactivate mold, bacteria, viruses, fungi and any other microorganisms from the irradiated surfaces,
- **Suitable for:**
  - Hospitals (surgery units, treatment rooms, emergency rooms, patients' wards, isolation rooms, A&E places, dirty linen rooms),
  - Clinics (doctor's and treatment rooms),
  - Laboratories,
  - Pharmacies,
  - The following industries: pharmaceutical, food-processing, cosmetics,
  - Hotel industry and gastronomy,
- **Advantages of NBVE 110/55 SL germicidal lamp:**
  - Mounting type: C - ceiling,
  - L - inductive counter mounted as standard on this model of germicidal lamp. The display indicates the time with an accuracy of up to 1 hour,
  - LW - electronic counter mounted on special order. Four-pole type display with indication to the hour. This model is also equipped with the function of acoustically informing the user of the need to replace the radiator,
  - Highest disinfection parameters confirmed by tests in microbiological laboratories,
  - The device has been awarded with numerous certificates,
  - High energy-saving properties,
  - Method of operation of external radiators: similar to the use of direct-acting germicidal lamps, no one is allowed to be inside the room during the operation of the lamp with external radiators running. They emit harmful radiation to the environment,
  - Application: in waiting rooms, emergency rooms, operating rooms, laboratories, corridors, nurseries and kindergartens, and in all places where there are many people and where strict sterility is required,
  - The acid-resistant sheet metal used in the case manufacturing process increases the disinfecting properties of the lamp,
- **Technical data:**
  - Ceiling-mounted lamp,
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 145 W,
  - Internal UV-C emitting element: 2 x 55W,
  - External UV-C emitting element: 1 x 55W,
  - UV-C irradiance at 1 m distance: 150  $\mu$ W/cm<sup>2</sup>,
  - Radiant life: 9000 h,
  - Fan capacity: 199 m<sup>3</sup>/h,
  - Disinfected volume: 45 - 90 m<sup>3</sup>,
  - Operating range of the lamp: 18 - 36 m<sup>2</sup>,
  - Anti-schock safety class: I,
  - Case type: IP 20,
  - Chamber's dimensions: 1125 x 285 x 130 mm,
  - Dimensions: 1190 x 400 x 130 mm,
  - Weight: 9,5 kg,
- It has a valid Technical Passport,
- Guarantee: 12 months,

