



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



NBV 15 S germicidal lamp (New)

- A brand new product,
- Made in Poland,
- NBV 15 S direct-acting germicidal lamp without a timer,
- This device emits UV-C radiation with a wavelength of 253.7 nm,
- The lamp irreversibly inactivates bacteria, viruses, molds, fungi and all other microorganisms,
- **The use of the device:**
 - Doctors' offices,
 - Treatment rooms,
 - Patient rooms,
 - Operating rooms,
- **Advantages of NBV 15 S germicidal lamp:**
 - Mounting type: C - ceiling,
 - Used to raise and maintain the level of microbiological purity in the room,
 - The lamps emit UV-C radiation with a wavelength of 253.7 nm,
 - They irreversibly inactivate bacteria, viruses, molds, fungi and all other microorganisms,
 - The germicidal lamp completely disinfects the air and surface,
 - The lamp reflector is made of high-quality reflective aluminum,
 - They are used in operating rooms, surgery rooms, bandaging rooms, patient rooms, clinics, laboratories, pharmacies, pharmaceutical, food and cosmetic industries,
 - Effectiveness of germicidal lamps approved by tests,
 - All required certificates,
- **Technical data:**
 - Ceiling lamp,
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 18 W,
 - UV-C emitting element: 15 W,
 - UV-C irradiance at 1 m distance: 0.9 W / m²,
 - Radiant life: 9000 h,
 - Disinfectable area: 6 - 8 m²,
 - Lamp rotation (exposure angle can be adjusted): 230°,
 - Anti-schock safety class: I,
 - Case type: IP 20,
 - Operation type: continuous,
 - Chamber's dimensions: 500 x 85 x 135 mm,
 - Chamber's weight: 2 kg,
 - The extension's length: 500 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

