



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **NBV 15 S germicidal lamp (New)**

- A brand new product,
- Made in Poland,
- NBV 15 S direct-acting germicidal lamp without a timer,
- This device emits UV-C radiation with a wavelength of 253.7 nm,
- The lamp irreversibly inactivates bacteria, viruses, molds, fungi and all other microorganisms,
- **The use of the device:**
  - Doctors' offices,
  - Treatment rooms,
  - Patient rooms,
  - Operating rooms,
- **Advantages of NBV 15 S germicidal lamp:**
  - Mounting type: C - ceiling,
  - Used to raise and maintain the level of microbiological purity in the room,
  - The lamps emit UV-C radiation with a wavelength of 253.7 nm,
  - They irreversibly inactivate bacteria, viruses, molds, fungi and all other microorganisms,
  - The germicidal lamp completely disinfects the air and surface,
  - The lamp reflector is made of high-quality reflective aluminum,
  - They are used in operating rooms, surgery rooms, bandaging rooms, patient rooms, clinics, laboratories, pharmacies, pharmaceutical, food and cosmetic industries,
  - Effectiveness of germicidal lamps approved by tests,
  - All required certificates,
- **Technical data:**
  - Ceiling lamp,
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 18 W,
  - UV-C emitting element: 15 W,
  - UV-C irradiance at 1 m distance: 0.9 W / m<sup>2</sup>,
  - Radiant life: 9000 h,
  - Disinfectable area: 6 - 8 m<sup>2</sup>,
  - Lamp rotation (exposure angle can be adjusted): 230°,
  - Anti-schock safety class: I,
  - Case type: IP 20,
  - Operation type: continuous,
  - Chamber's dimensions: 500 x 85 x 135 mm,
  - Chamber's weight: 2 kg,
  - The extension's length: 500 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

