



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **NBV 30 N germicidal lamp (New)**

- A brand new product,
- Made in Poland,
- NBV 30 N direct-acting germicidal lamp without a timer,
- This device emits UV-C radiation with a wavelength of 253.7 nm,
- The lamp irreversibly inactivates bacteria, viruses, molds, fungi and all other microorganisms,
- **The use of the device:**
  - Operating rooms,
  - Treatment and bandaging rooms,
  - Patient rooms,
  - Clinics,
  - Laboratories,
  - Pharmacies,
  - Pharmaceutical industry,
  - Food industry,
  - Cosmetics industry,
- **Advantages of NBV 30 N germicidal lamp:**
  - Type: N - universal, wall and ceiling direct acting germicidal lamp,
  - Purpose: to raise and maintain the level of microbiological cleanliness in rooms of various types,
  - Purpose: irreversible deactivation of all kinds of bacteria, viruses, molds, fungi and microorganisms,
  - Effect: effective disinfection of air and exposed surface,
  - Purpose: the use in operating rooms, treatment and bandaging rooms, operating rooms, clinics, laboratories and pharmacies, as well as in the pharmaceutical, food and cosmetic industries,
    - Additional information: this lamp is equipped with a high-quality aluminum reflector; in addition, it has all the required certificates, and its effectiveness is confirmed by tests,
- **Technical data:**
  - Wall mounted lamp,
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 33 W,
  - UV-C emitting element: 1 x 30 W,
  - UV-C irradiance at 1 m distance: 2.3 W / m<sup>2</sup>,
  - Radiant life: 9000 h,
  - Disinfectable area: 12 - 15 m<sup>2</sup>,
  - Lamp rotation (exposure angle can be adjusted): 200°,
  - Anti-schock safety class: I,
  - Case type: IP 20,
  - Operation type: continuous,
  - Chamber's dimensions: 960 x 85 x 135 mm,
  - Chamber's weight: 3 kg,
  - The extension's length: 120 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

