



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NBV 2 x 30 S germicidal lamp (New)**

- A brand new product,
- Made in Poland,
- Direct-acting germicidal lamp NBV 2x30 S: without timer,
- Irreversibly deactivates viruses and microorganisms,
- Equipped with a reflector made of reflective aluminum, which guarantees the highest quality of the device,
- **The use of the device:**
  - Operating rooms,
  - Treatment and bandaging rooms,
  - Patient rooms,
  - Clinics,
  - Laboratories,
  - Pharmacies,
  - Pharmaceutical industry,
  - Food industry,
  - Cosmetics industry,
- **Advantages of NBV 2x30 S germicidal lamp:**
  - Mounting type: C - ceiling,
  - These products are used to raise and maintain the level of microbiological cleanliness of interiors,
  - The lamps emit UV-C radiation with a wavelength of 253.7 nm,
  - They irreversibly inactivate bacteria, viruses, molds, fungi and all other microorganisms,
  - The germicidal lamp effectively disinfects the air and surface,
  - The reflector was manufactured from high-quality reflective aluminum,
  - These products are used in operating rooms, treatment rooms, bandaging rooms, patient rooms, clinics, laboratories, pharmacies, pharmaceutical, food and cosmetic industries,
  - The operational reliability of germicidal lamps is confirmed by tests,
  - All required certificates,
- **Technical data:**
  - Ceiling mounted lamp,
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 66 W,
  - UV-C emitting element: 2 x 30W,
  - UV-C irradiance at 1 m distance: 3.6 W / m<sup>2</sup>,
  - Radiant life: 9000 h,
  - Disinfectable area: 18 - 22 m<sup>2</sup>,
  - Lamp rotation (exposure angle can be adjusted): 230°,
  - Anti-schock safety class: I,
  - Case type: IP 20,
  - Operation type: continuous,
  - Chamber's dimensions: 960 x 85 x 145 mm,
  - Chamber's weight: 5 kg,
  - The extension's length: 500 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

