



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NBV 15 NL germicidal lamp (New)

- A brand new product,
- Made in Poland,
- Direct-acting germicidal lamp NBV 15 NL with timer,
- The device has a counter with a microprocessor, without a display, with an acoustic indication when the lamp is replaced,
- It has a CL-02 reader to read the current status of the counter,
- **The use of the device:**
 - Medical institutions,
 - Hospitals,
 - Pharmaceutical industry,
 - Food industry,
 - Cosmetics industry,
- **Advantages of NBV 15 NL germicidal lamp:**
 - L-counter with a microprocessor, without a display with an acoustic indication of the moment of replacing the radiator. This type of counter is mounted in direct-acting lamps,
 - Mounting method: W - wall-mounted,
 - With this germicidal lamp you will raise the level of microbiological purity,
 - Emitted radiation: UV-C with a wavelength of 253.7 nm,
 - When using this lamp, irreversible deactivation of viruses, molds, fungi and microorganisms occurs,
 - Purpose: effective disinfection of air and surfaces,
 - Reflector: this element is made of aluminum, which guarantees its highest quality,
 - Purpose: this germicidal lamp is used in medical institutes, hospitals and clinics, as well as in pharmaceutical, food and cosmetic industries,
 - The effectiveness of this lamp has been confirmed by tests; it is also proven by numerous certificates,
- **Technical data:**
 - Wall mounted lamp,
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 18 W,
 - UV-C emitting element: 1 x 15 W,
 - UV-C irradiance at 1 m distance: 0.9 W / m²,
 - Radiant life: 9000 h,
 - Disinfectable area: 6 - 8 m²,
 - Lamp rotation (exposure angle can be adjusted): 200°,
 - Anti-schock safety class: I,
 - Case type: IP 20,
 - Operation type: continuous,
 - Chamber's dimensions: 500 x 85 x 135 mm,
 - Chamber's weight: 2 kg,
 - The extension's length: 120 mm,
 - Counter type: without a display,
- It has a valid Technical Passport,
- Guarantee: 12 months,



