



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NBV 15 PL germicidal lamp (New)

- A brand new product,
- Made in Poland,
- Direct acting germicidal lamp NBV 15 PL with timer without display,
- The device has a counter with a microprocessor, without a display, with an acoustic indication when the lamp is replaced,
- It has a CL-02 reader to read the current status of the counter,
- **The use of the device:**
 - Medical institutions,
 - Hospitals,
 - Pharmaceutical industry,
 - Food industry,
 - Cosmetics industry,
- **Advantages of NBV 15 PL germicidal lamp:**
 - Counter type: L-counter digital with a microprocessor, without a display. This model audibly signals the moment when the radiant tube should be replaced. This type of counter performs its function perfectly in direct-acting germicidal lamps,
 - Mounting method: M - mobile on a tripod,
 - The device has passed tests that confirm its highest quality (it has the required certificates),
 - The lamp emits UV-C radiation with a wavelength of 253.7 nm,
 - Bacteria, molds, fungi and viruses exposed to this lamp are deactivated - they are no longer dangerous to people,
 - With the help of this lamp it is possible to make effective disinfection of air and surfaces,
 - The reflector of this lamp is made of reflective aluminum,
 - This device is designed for the use in operating rooms, treatment and bandaging rooms, patient rooms, clinics, laboratories, pharmaceutical, food and cosmetic industries,
- **Technical data:**
 - The lamp is placed on a mobile tripod,
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 18 W,
 - UV-C emitting element: 1 x 15 W,
 - UV-C irradiance at 1 m distance: 0.9 W / m²,
 - Radiant life: 9000 h,
 - Disinfectable area: 6 - 8 m²,
 - Lamp rotation (exposure angle can be adjusted): 270°,
 - Anti-schock safety class: I,
 - Case type: IP 20,
 - Operation type: continuous,
 - Chamber's dimensions: 500 x 85 x 135 mm,
 - Chamber's weight: 2 kg,
 - The extension's length: 1800 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

