



## **NBV 30 SL germicidal lamp (New)**

- A brand new product,
- Made in Poland,
- Direct-acting germicidal lamp NBV 30 SL with timer without display,
- The device has a counter with a microprocessor, without a display with an acoustic indication of the moment of replacing the radiator,
- **The use of the device:**
  - Medical institutions,
  - Hospitals,
  - Pharmaceutical industry,
  - Food industry,
  - Cosmetics industry,
- **Advantages of NBV 30 SL germicidal lamp:**
  - L - digital counter containing a micro-processor, without a display, with acoustic signaling of the moment of replacement, radiant, mounted in direct-acting lamps,
  - Lamp mounting method: C - ceiling mounted,
  - Lamps are used to raise and maintain the level of microbiological purity in the room,
  - Lamps emit UV-C radiation with a wavelength of 253.7 nm,
  - They irreversibly inactivate bacteria, viruses, molds, fungi and all other microorganisms,
  - The germicidal lamp completely disinfects the air and surface,
  - The lamp reflector is made of high-quality reflective aluminum,
  - The lamps are used in operating rooms, treatment rooms, bandaging rooms, patient rooms, clinics, laboratories, pharmacies, pharmaceutical, food and cosmetic industries,
    - Effectiveness of germicidal lamps approved by tests,
    - All required certifications,
    - It is recommended to turn on the lamps in rooms of standard height (up to 3m) for about 2-8 hours a day, and for 15 to 20 minutes before treatments,
- **Technical data:**
  - Ceiling mounted lamp,
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 33 W,
  - UV-C emitting element: 1 x 30 W,
  - UV-C irradiance at 1 m distance: 2.3 W / m<sup>2</sup>,
  - Radiant life: 9000 h,
  - Disinfectable area: 12 - 15 m<sup>2</sup>,
  - Lamp rotation (exposure angle can be adjusted): 230°,
  - Anti-schock safety class: I,
  - Case type: IP 20,
  - Operation type: continuous,
  - Chamber's dimensions: 960 x 85 x 135 mm,
  - Chamber's weight: 3 kg,
  - The extension's length: 500 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

