



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Germicidal Lamp NBV 2 x 30 S LW (New)

- New,
 - Made in Poland,
 - LW - digital meter with the mikro-processor, with display 4-field LED, with the acoustic signaling of the moment of the exchange the UV bulbs. The possibility of independent connecting to bought earlier direct and flow germicidal lamps,
 - The way of the assembly:
 - P - passing on the tripod,
 - they are designed to elevation and maintenance the level of microbiology cleanness of the rooms,
 - these devices emit the radiation UV-C about the length of the wave 253,7 nm,
 - the bacteria deactivate, viruses, moulds, mushrooms and all different micro-organisms,
 - Technical Data:
 - Power supply: 230 V 50 Hz,
 - Power requirement: 40 VA,
 - Element emits UV-C radiation: 30 W,
 - Radiation uv-c intensity at the distance of 1 m: 2,3 W/m²,
 - Bulb service lamp: 8000 h,
 - Effective area of the lamp: 12 - 15 m²,
 - Turn of the lamp (the possibility of the angle of exposing): 230°,
 - Class of protect. against electric shock: I,
 - Cover type: IP 20,
 - Kind of the work: continous,
 - Dome dimensions: 925 x 85 x 135 mm,
 - Dome mass: 3 kg,
 - The lenght of the arm: 500 mm,
 - Type of the counter: with 4 field display,
- Technical Passport,
- 12-month warranty

