



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FOTOVITA FV-10L lamp (Large) (New)

- A brand new product,
- Made in Poland,
- Phototherapy is much less expensive than drug therapy,
- Moreover, it can be used even for pregnant women and young breastfeeding mothers,
- The use of Fotovita lamps is not associated with any side effects,
- **Advantages of Fotovita lamps phototherapy:**
 - Quick and effective removal of the effects and symptoms of autumn and winter depression,
 - Eliminating the effects of sadness, anxiety, apathy and bad mood,
 - Effective fight against insomnia, period bulimia and stress that accompanies premenstrual tension,
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 250 W,
 - Method of powering fluorescent lamps: by high frequency (HF) inverters,
 - Illumination intensity: 10000 lux at a distance of 50 cm; 5000 lux at a distance of 75 cm,
 - Life of fluorescent lamps: 10,000 operating hours,
 - Anti-shock protection class: 1,
 - Protection class: I,
 - Dimensions [mm]: height: 570; width: 345; depth: 130 (without support), 280 (with support),
 - Weight: 5,0 kg,
- Provided with a valid Technical Passport,
- Guarantee 60 months,

