



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medical equipment table TALIE 3 (New)**

New,

Made in Poland,

It has a stable steel construction covered with powder paint,

As standard, the table is equipped with two metal shelves with a raised edge, a cable bracket and castors, 2 of which have a brake.

TALIE 3 model has adjustable shelves,

Additional options are available (at extra charge):

Additional metal shelf,

Additional wire support,

Power strip (4 sockets),

Technical specifications:

Width: 50 cm,

Height: 92 cm,

Depth: 40 cm,

Shelf dimensions (W/D): 41 x 40 cm,

Dead weight: 22 kg,

Maximum load per shelf: 15 kg,

Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

