



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Ultrasound Scanner Eurosonic 15 VET (New)

- New,
- Made in Poland,
- Eurosonic is the add-on module that could be connected to any personal computer through the USB-2.0 interface. Full digital design is a key of an excellent image quality and the computer application allows to archive the most important examination data like images, examination description and cine-loop files on hard disk or any other electronic storage media,
- Scanner is equipped with the sector probe operating at 3,5/5/7,5MHz frequencies,
- Comes with:
  - Ultrasound Scanner module,
  - 3,5/5/7,5 MHz Probe,
  - USB cable,
  - CD with installation software and user manual, solid and light weight transportation case
- Technical Data:
- Probe: Multifrequency 3,5/5/7,5 MHz
- Presentation modes: B/B+2B/Real-time zoom,
- Measurements: Distance, perimeter, area, volume,
- Calculation programmes: fetus age estimation for canine and feline based on HD and GS method,
- Grey scale: 256,
- Postprocessing: Interframe filtration, bilinear filtration,
- Computer communication: USB 2.0,
- Penetration depth: 60/80/120/160 mm,
- Cine loop: 256 frames (~20 seconds of examination),
- Data archive: Images, cine-loop files including description, time and date of an examination on hard disk or external storage device,
- Power consumption: 5 VA,
- Dimensions: 180 x 127 x 41 mm,
- Case dimensions: 390 x 314 x 102 mm,
- Weight: 430 g,
- Weight including case: 2300 g,
- Zasilanie: USB interface,
- 24-month warranty

