



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### Ambulance Holder for LIFEPAK 12 DLOUHY (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Austria,
- Intended for LIFEPAK 12,
- Made from powder coated steel,
- Equipped with solid metal handles,
- Foldable,
- Cummed base,
- Dimensions: 40 x 28 x 38 (51) cm,
- Technical Passport,
- Warranty: 30-days (6-month prolongation is extra paid),
- **The price does not include the defib,**
- **Compatible with LIFEPAK 15 and LIFEPAK 12**

