



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Examination Chair AKRUS SC 5010 HS (New)**

New,  
German production,  
Power supply 230V,  
Frequency 50/60Hz,  
Surgical table designed for e.g:  
ENT surgery,  
Maxillofacial surgery  
Plastic surgery,  
Electrically adjustable headrest, backrest, seat and leg rest,  
Each element individually adjustable,  
Adjustable seat height - possibility of extremely low adjustment 64-100 cm,  
Seat level adjustment,  
Possibility to mount various ergonomic headrests,  
Mobile - on 125mm castors with brakes on both sides,  
Possibility of battery operation,  
Large selection of additional accessories,  
Dimensions:  
Length of wheel base 760 mm,  
Width of the caster 580 mm,  
Seat width plus armrest 700 mm,  
Seat width 590 mm,  
Backrest height with upright backrest 1700 mm,  
Max. Chair length, backrest horizontal 1750 mm,  
Max. Patient weight 200/300 kg,  
Chair weight (depending on options) 114 kg,  
Optional equipment (extra charge):  
Standard headrest,  
Multi articulating headrest,  
Flat headrest,  
Articulated arm for additional headrest +30/-50°,  
Screen hanger,  
Ventilated screen holder,  
Armrest, adjustable in height, pivoting, on ball head,  
Arm rest, not adjustable,  
IV-Pole,  
Side rail,  
Foot switch,  
Foot joystick 5 segments,  
Battery,  
Holds the latest Technical Passport,  
Warranty 12 months,  
Finance options (Poland only): Instalments, Leasing, Loan

