



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mobile Examination Chair AKRUS SC 5010 SEK (New)

New,
German production,
Power supply 230V,
Frequency 50/60Hz,
Mobile operating table designed for, among others:
Orthopaedics,
General surgery
Dermatology
Plastic surgery,
Detachable arm supports,
Two independent leg supports,
Electrical height adjustment from 680 to 980mm,
Backrest manually foldable,
Electric seat adjustment from 0° to 20°,
Possibility to work on a battery,
Mobile - on 125mm castors with brakes on both sides,
Transparent elements in AKRUS table model: SC 5010 SEK:
Split foot segment,
Table for hand operation,
Shoulder segment,
Dimensions:
Length of wheel base (backrest) 760 mm,
Width of the caster base 580 mm,
Width 700 mm,
Seat width 590 mm,
Height of the backrest with straight back 1440 mm,
Max. Length, backrest horizontal 1960 mm,
Max. Patient weight 135/250 kg,
Weight (depending on options) 95 kg,
Optional equipment (at extra charge):
Standard headrest,
Multi articulating headrest,
Double head stabilizer,
Sliding headrest, flat,
Screen hanger,
Ventilated screen hanger,
Knee clamp,
Electric knee clamp,
Goepel leg rests (pair),
Arm support for side or stomach position,
Armrest,
Arm support for elbow surgery,
Arm support for hand surgery,
Back or breast stabiliser for lateral position,
Drip stand,
Side rail,
Foot switch,
Rechargeable battery,
Holds a valid Technical Passport,
Warranty: 12 months,
Finance options (Poland only): Installments, Leasing, Loan

