



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Laser AMS LASERSCOPE GREEN LIGHT HPS (Reconditioned)

Used, perfect technical condition,
Production of U.S.A (2009),
Power supply: 240 V, (powered by one phase power socket),
Frequency: 50-60 Hz,
Power: 20 A,
GreenLightHPS Green Laser - It has opened a new era of healing millions of men around the world. In the western countries, the procedure of evaporation of the hypertrophied prostate gland, due to its effectiveness and safety, has become the gold standard for urological treatment. In Poland, this modern surgical method has already helped several thousand patients, and scientific observations confirm its effectiveness.
Laser type: fixed, dual frequency,
Wavelength: 532 nm,
Average power levels: 20-120W in increments of 10 W,
Target beam: laser diode, red, 635 nm,
Output beam deviation: 0.05-0.15 radian of full angle at half power,
Built-in display: Touch screen,
Retractable fiber optic support,
Included:
Central control unit GreenLight HPS laser,
Foot double attachment,
3 pieces of protective glasses,
2 x key for switching on the laser,
4 x LASEROSCOPEGreenLightHPS optical fiber:
Fiber GL HPS,
After the deadline,
Sterilization: single use,
Weight: 152 kg,
Dimensions: 42 x 86 x 110 cm,
Mobile on wheels with brakes,
User manual in Polish, German,
Has a current Technical Passport issued,
1-month warranty

