



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical laser OPTOTEK OPTOLIGHT 25 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Slovenia,
- Year of Manufacture: 2009,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 650 W,
- Laser power: 25 W,
- **OPTOLIGHT 25** - is a medical laser for soft tissue surgeries. Supplies up to 25 W through a 200µm fiber optics. It is perfect for high-precision surgery, such as: soft tissue surgery in dentistry, endoscope treatment or venous surgery.
- **Laser is used in:**
 - Dental treatment on soft tissue, f.e. perodontium treat with curretage method, gums modelling, frenum cut, buccal cavity anaplasty, cutting and removing tissue, sterilization at root canal treatment.
 - Endoscope treatment: cutting, coagulation or heat therapy of soft tissues.
 - Venous surgery for clogged veins.
- Built-in LCD display (128 x 240 p),
- **Technical data:**
 - Laser type: laser diode LASER class 4,
 - Radiation wavelength: 980 nm,
 - Laser power: 1-25 W, increment: 1W,
 - Work mode: continous, single impulse, repeated impulse,
 - Single impulse mode: time set between 25 ms - 30 s,
 - Repeated impulse mode:
 - Impulse time: 5 ms - 9 s.,
 - Pause time: 20 ms - 9 s.,
 - Lasera: laser sight 635 nm, power 1 mW, regulated intensity,
 - Transferring beam system: fibers 200 / 400 / 600 um / N.A 0,22-SMA standard connector 905 to a laser socket,
 - Foiled control panel,
 - Weight: 14 kg,
- **Set includes:**
 - Main device - laser,
 - Fiber optics in PVC cover:
 - Super Flex 200,
 - Sure Flex 365,
 - Sure Flex 150,
 - Foot pedal,
 - Protective glasses - 4 pairs,
 - Ignition key,
- User manual,
- Size: 300 x 395 x 250 mm
- Technical passport,
- Warranty: 30-days

