



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Pulmonary fan DRAGER RESPICARE VV (Reconditioned)

Used, very good technical condition,
French production,
Power supply: 230 V,
DC: 12 to 36 V,
Frequency: 50-60 Hz,
Power: 70 VA,
Easy to use; simple operation, low noise,
Only 5 buttons to adjust all parameters,
Liquid crystal display (with standby mode for energy saving), in which all parameters are visible, including curves and barcodes,
Keyboard lock to prevent accidental changes,
Works in any position,
Ventilation modes: IPPV (S) IPPV (P) SIMV Ech plateau pressure,
Current volume: from 0.1 to 1.6 liters,
Minute volume: 3 to 25 l / min,
Frequency: from 8 to 50 cycles / min; in IPPV (P) mode, 4 to 20 cycles / min.
I / E ratio: 1/1 to 1/3,
Trigger: - 1 to - 3 cm H₂O,
Insp. In PVC mode: 0.4 to 1.5 L / s,
Ech: available in IPPV and IPPV (S) mode,
Low pressure alarm: 2 to 40 cmH₂O,
High pressure alarm: 10 to 60 cm H₂O,
Software update via the RS 232 port,
Included:
Central unit DRAGER RESPICARE VV,
The hose with the mask,
Nasal cannula,
3 spare filters,
Overall dimensions: 240 x 255 x 385 mm,
Weight: about 12 kg,
Has a current Technical Passport issued.
1-month warranty

