



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Therapy Batriatric 2 Section Couches (ST2573-SEERSMEDICAL) rehabilitation table (New)

New,
English production,
The height of the table and the backrest is adjusted by means of an electric actuator,
Retractable wheels and adjustable feet provide excellent stability (RWD),
Frame construction designed for easy use of a mobile lift,
The table as a standard has a hand-held remote control,
Cylinders used in the TiMotion table,
Technical data:
Total length of the table bed: 200 cm,
The width of the bed is 80 cm,
Safe working load and lifting capacity (from minimum height) 325 kg (715 lbs)
Adjustment of the back section in the range from 0 ° to 75 °,
Table height adjustment range from 45 cm to 98 cm,
Handheld remote control and face hole with a standard plug,
Exclusive 80 cm upholstery as standard,
Additional accessories (extra charge):
Double foot switch (ST6207)
Has a current Technical Passport issued,
Warranty 24 - months for electric actuators



TAPICERKI STANDARD + [ZA DOPŁATĄ]



TAPICERKI STANDARD [W CENIE]

