



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Welding machine TECNO SEAL PRINT (New)**

New

Italian production,

Tecno Seal is a roller sealer with an automatic packaging system,

The welding machine uses the latest technology used for medical seals, used in medical centers, dental offices, hospitals and clinics.

Advantages of the welding machine:

Guarantees very high output power, maximum operating speed without routine maintenance; applies bags with 3 sides sealing (available in various sizes,

Another feature is an electronic control card - enabling the management of usage parameters,

The use of this device is simple and rational - the operator only has to put a bag containing the device to support the welding machine;

photocell automatically, autonomously feeds the bag and seals it,

Integrated printer included (Art. TA002TAU)

Technical data:

Overall dimensions:

Width: 555 mm,

Height: 250 mm,

Depth: 145 mm,

Weight: 14 kg,

Power consumption 500 W,

Power supply: 230 V / 50/60 Hz,

The width of the seal: 12 mm,

Sealing length: unlimited,

Welding speed: 10 m / min,

Welding temperature: max. 220 ° C

Has a current Technical Passport issued,

12 months warranty

